

PC recommendation & specifications

At Nanolive we work with the Acer Predator, which gives very good user experience and performance. Therefore, for optimal performance, we recommend this specific computer to our clients.

You are free to choose the brand that you like, please be aware of the recommended technical specifications mentioned below to guarantee a smooth operation with our software STEVE.

ACER ASPIRE PREDATOR G3-605



The recommended technical specifications of this ACER PREDATOR are the following:

Operating System: 64-bit versions of Windows Vista, 7 or 8
CPU: Intel core i7 3GHz
Memory: 8GB DDR3 RAM
Hard Drive: 1 TB HDD
GPU: NVIDIA GeForce GTX 970 4GB

(minimum requirement spending evaluation)

Enjoy exploring living cells now!

Your Nanolive Team

Nanolive SA